



Model	Sand Rider
Overall Length	72"
Overall Width	38"
Overall Height	48"
Assembled weight	49 pounds
Front wheel	11.8"
Rear wheels	19.3"
Recommended Maximum Weight	350 pounds
Seat height	19"
Seat width	18"
Seat depth	16"
Inter armrest spacing (standard)	16"
Inter armrest spacing (comfort fit)	21"

## Materials:

**Frame:** Powder coated\* aluminum 1.7" diameter frame, 1.3" diameter push handle and armrests

**Upholstery\*:** Sling quality vinyl coated polyester, UV stable, Microban® antimicrobial protection, GREENGUARD® Select Certified

**Fasteners:** 316 stainless steel bolts and nuts with nylon washers

**Wheels:** Wheeliez™ low pressure polyurethane tires, polypropylene hub, 1" bearing

(\*Custom units- upholstery fabrics and frame coating materials may vary on request)